

THE LABORATORY ANIMAL PORT AUTHORITY

How Low Can You Go



ULTRA-LOW DEAD VOLUME IMPLANTED PORTS FOR STUDIES INVOLVING CSF AND LOW-VOLUME TEST ARTICLES

The LoVol™ ports are another innovation from Instech Solomon. The ports comprise ultra-low dead volumes to enhance CSF access and delivery of costly low-volume test articles. Use in conjunction with a smaller catheter to achieve the smallest dead volumes.



SPECIFICATIONS	MINLO	MIDLO
Dead Volume	15µL	40µL
Body Material	Titanium	Titanium
Height	.275 in (7.0 mm)	.395 in (10.0 mm)
Weight	3.0g	6.9g
Catherter Sizes	3Fr-5Fr	3Fr-5Fr

The ultra-low dead volume is packaged into a conventional SoloPort™ body. These LoVol™ products are based on the SoloPort™ MIN and MID port designs.



As such, your accessing techniques are unchanged, including use of 24g or 22g Huber needles.

**INSTECH
SOLOMON**

www.instechlabs.com

5209 Militia Hill Road
Plymouth Meeting, PA 19462
USA

© 2011 Solomon Scientific

TL + 1 800 443 4227
TL + 1 610 941 0132
FX + 1 610 941 0134

All rights reserved

SOLOMON
S C I E N T I F I C

www.solsci.com