

Product

MODEL 100 SYRINGE PUMP

Flow rate 0.01 – 3,000.00 mL/hr

Accuracy \pm 3% (see B-D statement on syringe accuracy)

Dimensions 9.5 x 4.0 x 5.9 inches / 24 x 10 x 15 cm

Weight 4.9 pounds (2.2 kg)

AC power Order either 110V (US cord) or 220V (Europe cord)

DC Battery backup 12 hours approx.; <8 hour recharge time

Syringe size 1 cc – 140 cc plastic and glass

Syringe brands B-D® (Becton Dickinson), Monoject® (Covidien), Terumo

Regulatory 510(k) (human-use), CE Marking (non-human use)

Made in US

Compatible with Orchestra™ software and wireless network